

MODEL

| | |
|---|---|
| Jazzy Elite ES-2S-SS ¹ | Group 2: Standard (K0822) |
| Jazzy Elite ES-2S-C ¹ | Group 2: Standard (K0823) |
| Weight Capacity | 300 lbs. |
| Turning Radius | 24.75" |
| Overall Width³ | 22.75" |
| Overall Length³ | 34.75"/41" (with footplate) |
| Maximum Speed² | Up to 4 mph |
| Ground Clearance³ | 1.85" at motors |
| Front Wheels | 3" solid anti-tips |
| Drive Wheels | 9" solid |
| Rear Wheels | 6" solid casters |
| Drivetrain | Two-motor, in-line, front-wheel drive |
| Braking System | Regenerative and electro-mechanical |
| Suspension Type | Limited |
| High-back Seat Sizes | 16' x 16"-18", 18" x 18"-20", 20" x 18"-20" |
| Seat-to-Floor Range | 21.5" - 23.5" |
| Per Charge Range^{3,5} | Up to 15 miles |
| Standard Electronics | 40A, PG GC3 |
| Component Weights | Base 69 lbs. (<i>heaviest piece</i>) Seat 50.2 lbs. (<i>comfort high-back</i>) Batteries 25 lbs. each |
| Battery Requirements^{4,6} | (2) 12 volt, U-1 batteries |
| Battery Charger | 3.5 amp, off-board |
| Battery HCPCS | E2365 |
| Warranty | Frame: 5-year limited Drivetrain: 13-month Electronics: 13-month Battery: 13-month |

¹ Jazzy® Power Chairs are FDA Class II Medical Devices designed to aid individuals with mobility impairments.

² Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).

³ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

⁴ AGM or gel cell type required.

⁵ Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

⁶ Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

(US) 800-800-8586 • (Canada) 888-570-1113

pridemobility.com |    

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

The information contained herein is correct at time of publication; we reserve the right to alter specifications without prior notice.

0004461-52322

Jazzy^{® 1}
ELITE ES

**OPTIMAL INDOOR
MANEUVERABILITY
IN SMALL SPACES**



FDA Class II
Medical Device¹

AMERICA'S FAVORITE POWER CHAIR

pridemobility.com |    

POWER, PERFORMANCE AND STYLE



Depth-adjustable, limited-recline, fold-flat, high-back seat

Height, width and angle adjustable armrests

Intuitive controller is easy to operate

Flip-up foot platform offers greater comfort

FDA Class II Medical Device¹

An experience so superior, it can only be elite. The **Jazzy® Elite ES¹** features large, front-drive wheels to absorb the bumps and uneven surfaces of everyday life. The comfortable, high-back seat gives you the support you need all day, every day.

- Depth-adjustable, limited-recline, fold-flat, high-back seat
- Universal, fixed, in-line joystick mount
- Flip-up foot platform
- Black, flat-free, non-marking tires
- Easy-access free-wheel lever
- 60" lap belt

Colors



Candy Apple Red

Seating

Black, vinyl, depth-adjustable and limited recline standard comfort high-back with 60" lap belt available sizes:

- 16"W x 16" - 18"D
- 18"W x 18" - 20"D
- 20"W x 18" - 20"D

Options

- Bell-shaped foot platform
- Elevating leg rests
- Heel loops
- Height-adjustable foot platform bracket
- Residual limb support
- Swing-away legrests
- Swing-away joystick

Optional Seating

- Black vinyl, high-back seat up to 20" x 20"

Accessories

- Cane/crutch holder
- Cell phone holder
- Cup holder
- Flag holder
- Lap belt (50", 70", 80", 90")
- Oxygen tank holder
- Rear basket
- Saddle bag
- Walker holder
- Weather cover
- XLR USB charger



Front-wheel drive design



Easy-access freewheel levers



PG GC3 Controller

Military Patches

