# 1 FOR REHAB POWER



# Q6 edg<sub>H</sub>

With features like standard 4-pole motors, Mid-Wheel 6® Drive Design and

ATX Suspension, the Q6 Edge® HD is engineered to meet the performance needs of the most active user.

The Q6 Edge HD accepts a complete range of seating and electronics options to deliver exceptional rehab capability at an excellent value, and is available with our iLevel® technology.



## Specifications

<u> </u>	
DRIVE WHEELS	Standard: 14" Run flat Optional: 14" Pneumatic 14" Knobby
CASTER WHEELS	Front: 6" Solid Rear: 6" Solid
ANTI-TIP WHEELS	N/A
MAXIMUM SPEED <sup>2</sup>	Up to 4.7 mph
GROUND CLEARANCE	3" (frame center)
TURNING RADIUS 3	19.75"
OVERALL LENGTH	35.6" without foot riggings
BASE WIDTH	26.25"
SEATING SIZES 16	
Synergy®/Static Seat	Width: 20-32" Depth: 20-28"
TRU-Balance® 3 HD Power Tilt	Width: 20-32" Depth: 20-25"
TRU-Balance® 3 HD Lift & Tilt	Width: 20-24" Depth: 20-25"
SEAT-TO-FLOOR HEIGHTS	
Synergy®/Static Seat <sup>6</sup>	16.5-18.5"
TRU-Balance® 3 HD Power Tilt <sup>6</sup>	18.5-20.5"
TRU-Balance® 3 HD Lift & Tilt <sup>6</sup>	17.75-19.75"
MANUAL TILT/RECLINE <sup>7</sup>	No
BATTERY SIZE <sup>4</sup>	Group 24
BATTERY WEIGHT <sup>9</sup>	52 lbs.
AVAILABLE ELECTRONICS 8	120A Q-Logic 3 EX 120A Q-Logic 3e 75A Q-Logic NE+ 75A Q-Logic NE
BATTERY CHARGER	8A Off-board
MOTOR PACKAGES	4-Pole 4.7 mph <sup>1</sup>
WEIGHT CAPACITY 5	450 lbs. <sup>10</sup>
BASE WEIGHT	143.5 lbs.
MEDICARE CODE COVERAGE	<b>Group 3</b> : HD (K0850/K0851), Single Power HD (K0858), Multiple Power HD (K0862)

- 1) Standard
- 2) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 3) Without front riggings
- 4) Battery capacity may change when adding power positioning systems.
- 5) Weight capacity is dictated by the selection of seating.
- 6) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless
- The Synergy® Manual Tilt weight capacity is 75 lbs. and maximum dimensions are 14" W x 14" D. Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details. All Quantum Rehab® units can be ordered with specialty controls
- Battery weight may vary +/- 3 lbs. based on manufacturer.
- 10) Weight capacity is reduced to 400 lbs. with the iLevel® power adjustable seat height
- 11) A minimum seating size of 20"x20" is required with iLevel®. iLevel® maximum seating size is 24" x 25".

#### **Features**

- Mid-Wheel 6® allows six wheels on the ground for maximum stability
- Low impact OMNI-Casters (nylon, spherical-shaped casters) on front and rear prevent wheel hang-ups
- ATX Suspension (Active-Trac® with extra stability) incorporates front OMNI-Casters and semi-independent rear caster beam for enhanced performance over more varied terrain
- iLevel® option enables up to 10" of seat elevation while at walking speed (up to 3.5 mph)<sup>10,11</sup>

## **Options**

- iLevel® power adjustable seat height<sup>10,11</sup>
- TRU-Balance® 3 (450 lbs.) HD Tilt, HD Recline, and HD Tilt & Recline
- Q-Logic 3 Drive Control System
- USB mobile device charger
- Swing-away joystick mount
- LED fender lights (standard with iLevel®)

### Colors

SEAT BACKS for units without TRU-Balance® 3 are only available in Back in Black.



Pretty



in Pink

Green Machine

Purple

Evolution Orange







Electric Blue

Red Corvette

Grey Street

Midnight Blue





Night Train

Back in Black





