



Revolving Around Your Needs.



Wheelchairs
Complete Selection. Clinical Expertise.

High Quality, Ready To Roll.

The durability of our wheelchairs is famous, and we stand behind them with long warranties: one year on parts and lifetime on the frame. TIG welding strengthens the joints, an extra layer of upholstery support prevents sagging, brass inserts reduce stripping from upholstery changes, and a coating on the frame helps you prevent rust and maintain a brand-new appearance.



MDS806250EV

Excel K2B Basic Wheelchairs

- » Carbon steel frame with rugged silver-vein finish
- » Upholstered, padded armrests
- » Dual-axle hemi-height adjustable
- » Chart pocket on the back
- » Easy-to-clean vinyl upholstery
- » Desk-length removable arms
- » 300-lb. (136 kg) weight capacity
- » Warranty: 12 months on parts, lifetime on the frame and welding

| Item No. | Leg Rests/Footrests | Seat Width | Color |
|--------------|----------------------|---------------|-------|
| MDS806250NEV | Swing-Away Footrests | 16" (40.6 cm) | Black |
| MDS806300NEV | Elevating Leg Rests | 16" (40.6 cm) | Black |
| MDS806250EV | Swing-Away Footrests | 18" (45.7 cm) | Black |
| MDS806300EV | Elevating Leg Rests | 18" (45.7 cm) | Black |
| MS806400EV | Swing-Away Footrests | 20" (50.8 cm) | Black |
| MDS806450EV | Elevating Leg Rests | 20" (50.8 cm) | Black |



MDS806300D

Excel 2000 & Excel 2000 Narrow Wheelchairs

- » Carbon steel frame with rust- and chip-resistant chrome plating
- » Upholstered, padded armrests
- » Dual-axle hemi-height adjustable
- » Chart pocket on the back
- » Easy-to-clean vinyl upholstery
- » Desk-length removable arms
- » 300-lb. (136 kg) weight capacity
- » Warranty: 12 months on parts, lifetime on frame and welding

| Item No. | Leg Rests/Footrests | Seat Width | Color |
|------------|----------------------|---------------|-------|
| MDS806250N | Swing-Away Footrests | 16" (40.6 cm) | Navy |
| MDS806300N | Elevating Leg Rests | 16" (40.6 cm) | Navy |
| MDS806250D | Swing-Away Footrests | 18" (45.7 cm) | Black |
| MDS806300D | Elevating Leg Rests | 18" (45.7 cm) | Black |
| MDS806400 | Swing-Away Footrests | 20" (50.8 cm) | Black |
| MSD806450 | Elevating Leg Rests | 20" (50.8 cm) | Black |



MDS806850

Extra-Wide Wheelchairs

- » Carbon steel frame with rust- and chip-resistant chrome plating
- » Upholstered, padded armrests and calf pads
- » Chart pocket on the back
- » Easy-to-clean vinyl upholstery
- » Desk-length removable arms
- » 500-lb. (226.8 kg) weight capacity
- » Warranty: 12 months on parts, lifetime on the frame and welding

| Item No. | Leg Rests/Footrests | Seat Width | Color |
|-----------|----------------------|---------------|-------|
| MDS806700 | Swing-Away Footrests | 20" (50.8 cm) | Navy |
| MDS806750 | Elevating Leg Rests | 20" (50.8 cm) | Navy |
| MDS806800 | Swing-Away Footrests | 22" (55.9 cm) | Navy |
| MDS806850 | Elevating Leg Rests | 22" (55.9 cm) | Navy |
| MDS806900 | Swing-Away Footrests | 24" (61 cm) | Navy |
| MDS806950 | Elevating Leg Rests | 24" (61 cm) | Navy |



MDS808450

Infinite Position Reclining Wheelchairs

- » Infinite position reclining goes from 90° to 140° with quick adjust handles
- » Full back and head support is provided by 22" (56 cm) high back and 10" (25.4 cm) removable headrest
- » Carbon steel frame with chip-resistant chrome plating
- » Durable vinyl upholstery
- » Hemi-height adjustable axles raise/lower seat 2" (5.08 cm)
- » Desk length removable arms
- » Anti-tippers standard
- » 300-lb. (136 kg) and 350-lb. (159 kg) weight capacities
- » Warranty: 12 months on parts, lifetime on frame and welding

| Item No. | Leg Rests/Footrests | Seat Width | Weight Capacity |
|-----------|---------------------|-------------|------------------|
| MDS808350 | Elevating Leg Rests | 16" (41 cm) | 300-lb. (136 kg) |
| MDS808450 | Elevating Leg Rests | 18" (46 cm) | 300-lb. (136 kg) |
| MDS808550 | Elevating Leg Rests | 20" (51 cm) | 350-lb. (159 kg) |
| MDS808650 | Elevating Leg Rests | 22" (56 cm) | 350-lb. (159 kg) |

Our Warranty

We stand behind all our wheelchairs with the industry's best warranty: 12 months on parts and lifetime on the frame and welding.



MDS809750

Shuttle Bariatric Wheelchair

- » Our strongest and widest wheelchair
- » Double cross-brace for maximum support and durability
- » Arms and legrests are removable to make exiting easy
- » Comfortable nylon upholstery and padded legrests
- » Desk length arms
- » Hemi-height adjustable with a triple-axle design
- » 850-lb. (385.6 kg) weight capacity
- » Warranty: 12 months on parts, lifetime on the frame and welding

| Item No. | Leg Rests/Footrests | Seat Width |
|-----------|---------------------|---------------|
| MDS80650 | Elevating Leg Rests | 24" (61 cm) |
| MDS809750 | Elevating Leg Rests | 26" (66 cm) |
| MDS809850 | Elevating Leg Rests | 28" (71.1 cm) |
| MDS809950 | Elevating Leg Rests | 30" (76.2 cm) |



MDS808200BAR

Freedom Plus Transport Wheelchair

- » Weighs just 33-lb. (15 kg)
- » Back folds down for easy storage and transport
- » Handbrakes and a seat belt for safety and control
- » Easy-to-clean vinyl upholstery
- » Anti-tippers
- » Push bar with foam grips, cup holder and a chart pocket on the back
- » Full-length permanent arms
- » 400-lb. (181 kg) weight capacity
- » Warranty: 12 months on parts, lifetime on the frame and welding

| Item No. | Leg Rests/Footrests | Seat Width |
|--------------|---------------------|---------------|
| MDS808200BAR | Elevating Leg Rests | 22" (55.9 cm) |

Our Warranty

We stand behind all our wheelchairs with the industry's best warranty: 12 months on parts and lifetime on the frame and welding.

Excel Kidz Pediatric Wheelchairs

- » Adjustable telescoping handles for caregiver comfort
- » Treaded, flat-free tires for safe riding
- » Anti-tippers included
- » Flip-back desk-length arms
- » Detachable elevating leg rests
- » Upholstered, padded armrests and calf pads
- » Dual-axle hemi-height adjustable
- » Easy-to-clean vinyl upholstery
- » 250-lb. (113 kg) weight capacity
- » Warranty: 12 months on parts, lifetime on the frame and welding



MDS806140PD

| Item No. | Leg Rests/Footrests | Seat Width |
|-------------|---------------------|---------------|
| MDS806140PD | Elevating Leg Rests | 14" (35.6 cm) |

Optional Wheelchair Accessories

| | | Excel K2B Basic | Excel 2000 | Extra-Wide | Reclining | Shuttle Bariatric | Freedom Plus Transport | Excel Kidz Pediatric |
|--|-------------|-----------------|------------|------------|-----------|-------------------|------------------------|----------------------|
| Description | Item No. | | | | | | | |
| Overhead Anti-Theft Device for 18" - 22" (45.7 - 55.9 cm) Chairs | MDS85197A | ● | ● | ● | | | | |
| Overhead Anti-Theft Device for 24" (61 cm) Chairs | MDS85197SH | | | ● | | ● | | |
| Anti-Fold Device for 16" - 20" (40.6 - 50.8 cm) Chairs | MDS85196 | ● | ● | ● | ● | | | |
| Anti-Fold Device for 22" - 24" (55.8 - 61 cm) Chairs | MDS851965SH | | ● | ● | ● | ● | ● | |
| Anti-Tip Device | MDS85189 | ● | ● | ● | | | | |
| Anti-Tip Device | MDS85189SH | | | | | ● | | |
| Anti-Tip Device | MDS85189FT | | | | | | ● | |
| O ₂ / IV Combo | MDS85190 | ● | ● | ● | | ● | | ● |
| O ₂ / IV Combo | MDS85190FT | | | | | | ● | |
| O ₂ Holder | MDS85181U | ● | ● | ● | | | | ● |
| O ₂ Holder | MDS85181FT | | | | | | ● | |
| IV Pole | MDS85183 | ● | ● | ● | | ● | | ● |
| IV Pole | MDS85183FT | | | | | | ● | |
| Tinerman Nuts | WCA806991 | ● | ● | ● | ● | ● | ● | |

Helping You Make Informed Purchasing Decisions

ONLINE TOOLS

Capital Quote Tool

www.medlinecapitalquote.com

It's an online experience that simplifies the way you search for and order equipment. You can also request a quote, and you'll receive a response within 24 hours from our team of specialists.

Wheelchair Selector

www.medline.com/pages/wheelchair-product-selector

Decide which wheelchair model is best for your facility with ease with our Wheelchair Selector Tool.

It asks you questions and makes suggestions of wheelchairs that fit your needs based on your answers.



Online Capital Quote Tool With Nurses Station

SUPPORT PROGRAMS

Wheelchair Assessments

At no cost to you, we can evaluate and analyze your wheelchair fleet. You receive a list of parts that you need for repairs and a list of wheelchairs that are past their useful lifespan, beyond repair or possibly dangerous. Contact your Medline representative to schedule an assessment.



Custom Laminated Label

Custom Labeling

With Medline's Custom Labeling program, you can advertise your facility's brand and specify department assignments, while deterring pilferage. Custom labels or stencils are free with your purchase of 10 or more wheelchairs. Custom embroideries are free with your purchase of 50 or more wheelchairs. For more information, contact your Medline representative.

Contact your Medline representative or call 1-800 MEDLINE to order yours today.



Medline Industries, Inc. One Medline Place, Mundelein, IL 60060

Medline United States
1-800-MEDLINE (633-5463)
medline.com | info@medline.com

Medline Canada
1-800-396-6996
medline.ca | canada@medline.com

Medline México
01-800-831-0898
medlinemexico.com | mexico@medline.com

FOLLOW US [f](#) [in](#) [t](#) [BLOG](#)