

LECKEY[®]

Distributed by  Mobility

mygo Stander



Encouraging a natural standing posture with improved function is as easy as **1-2-3**

Bringing standing therapy to children ages 4 to 14 years old as they encounter growing postural changes. As a child grows out of the popular Squiggles Stander, their standing therapy can continue with the **Mygo Stander** as it encourages a natural standing posture with improved function.

1 Stander

Mygo Stander offers continuous angle adjustments from vertical to 10 degrees from horizontal. Through the adjustment of the knee and foot supports, **Mygo Stander** can also accommodate up to 25 degree fixed hip and knee contractures.



2 Sizes

Available in two sizes to accommodate ages 4-14.

Size 2
8-14 years

Size 1
4-10 years



3 Configurations

A 3-in-1 stander offering prone, supine and upright standing in one product. The **Mygo Stander** provides clinicians with the tools to maximize posture and positioning.



Trunk Support

The wings on the contoured shoulder support can be angle adjusted to provide shoulder protraction. This brings the hands to midline which increases upper limb and hand control for improved function.

The ergonomically designed chest pad replicates the shape of the rib cage allowing the shoulder girdle to freely protract and retract in prone standing, improving upper body function.



Pelvic Stability

The pelvic derotation belt enables the therapist to position the client's pelvis in a central and neutral position with minimum effort.



Leg & Foot

The angle of the knee and foot supports relative to the hip and chest supports can be adjusted by 25° to accommodate contractures. The footplates can also be adjusted individually to accommodate up to 10° plantarflexion or dorsiflexion.



Specifications

	Size 1	Size 2
Age	4-10	8-14
User Height	42.5" - 57.1"	53.1" - 66.9"
Weight Capacity	110 lbs.	132 lbs.
Product Weight	48 lbs.	55 lbs.
Prone Range (0°=vertical)	0° - 80°	0° - 80°
Supine Range (0°= vertical)	0° - 80°	0° - 80°
Available Colors	Blue, Pink, Orange, Infection Control (Black)	

For complete specifications and available options, visit www.kimobility.com

LECKEY®