



K4 Basic Lightweight Wheelchair
K4 Basic Lightweight Wheelchair with Removable Swing-Back Desk-Length Arms and Swing-Away Footrests, 300 lb. Weight Capacity, 18" Width

- Durable, TIG-welded frame with gray powder coat finish plus swing-back, desk length arms with flip-back arm release for easy in and out
- Breathable nylon upholstery plus custom adjustability (dual axle hemi-height and back height) provides a more comfortable patient experience
- Smooth-rolling, solid, flat-free latex tires; quick release axles for easy, tool-free removal of rear wheels
- Optional accessories: O2 Holder (item MDS85181U), Anti-Fold Device (item MDS85196), IV Pole (MDS85183), Overhead Anti-Theft Device (MDS85197A), O2/IV Combo (MDS85190), Anti-Tip Device (MDS85189WK4), Tinnerman Leg Rest Locks (WCA806991)

Packaging	1 Case / Each
Manufacturer	Medline
Manufacturer #	MDS806500E

[Documents/SDS](#)

No SDS content is available for this product.

SKU MDS806500E

Specifications

ARM LENGTH	Desk Length
ASSEMBLY REQUIRED	No
AXLE TYPE	Dual
BACK HEIGHT	17 "
BARIATRIC	No
CASTER SIZE	8 "
COLOR	Black
FABRIC TYPE	Nylon
FOLDABLE	Yes
HPIS CODE	630_20_50_20
LATEX FREE	No
LEG REST STYLE	Swing Away
MRI SAFE	No
MATERIAL	Steel
OVERALL HEIGHT	38 "
OVERALL WIDTH	25.5 "
OVERALL WIDTH INCHES	25.5"
PRODUCT TYPE	Lightweight Wheelchair
PRODUCT WEIGHT (LB)	33.5 lb
QUICK RELEASE AXLES	No
REAR WHEEL SIZE	24 "
REMOVABLE ARM	Yes
SEAT DEPTH	16 "

SEAT WIDTH	18 "
SEAT TO FLOOR HEIGHT	20"
SELF PROPELLING	Yes
UNSPSC	42192210
WEIGHT CAPACITY	300 lb.

Unit of Measure	Conversion	Net/Gross Weight (lbs)	Volume (cubic ft)	Shipping Dimensions (inch) L x W x H	GTIN
Case (CS)	1.0 Each	0.0 / 45.415	8.825	34.5 x 13.0 x 34.0	80196305738
Each (EA)	1.0 Each	0.0 / 43.629	9.434	31.89 x 11.811 x 32.283	80196305738

California Customers Only... WARNING

HCPCS Information

HCPCS Code: K0004

HCPCS Disclaimer: HCPCS codes and Home Health Consolidated Billing codes provided by Medline are intended as general guidelines only. Medline does not guarantee coverage or reimbursement of any products. You must address all coverage and reimbursement issues (including the correctness and accuracy of codes) with your individual payers. It is your responsibility to ensure the accuracy and appropriateness of each claim you submit, in accordance with all applicable payer requirements.

For more information on HCPCS codes, Click here