

FRONT WHEEL DRIVE FOR GREAT MANEUVERABILITY AND TERRAIN HANDLING PLUS TOP PERFORMANCE AND STYLE









FEATURES

- Standard comfortable Captain's Seats with deep-contoured cushions and backs for maximum comfort, superior user positioning and lateral support.
- Captain's Seats swivel and feature a fold-down backrest.
- Adjustable height headrest.
- Seats are height adjustable, from 20.5" to 23".
- Padded, flip-up armrests allow for easy transfers and user access.
- Armrests are height and width adjustable for positioning, support and adaptability to the needs of the user.
- Adjustable length controller mount can be affixed to either the left or the right armrest and can be fitted with an optional swing-away arm for additional flexibility.
- Standard footplate can be flipped up out of the user's way and adjusted in both height and angle.
- Chair can be fitted with optional elevating legrests allowing the user to travel comfortably with rasied legs-a necessity for those with some medical conditions such as edema.
- Programmable controller can be personalized for individuals with special requirements.
- Flat-free, non-marking tires offer worry-free travel.
- Motor assemblies offer easy freewheel operation, allowing the chair to be pushed manually, which can be helpful for transport.
- Front anti-tip wheels for added safety.

COLORS

All Titan power chairs include interchangeable color panels in Blue and Red.



SEATING OPTIONS

18" (w) x 18" (d) Low Back Captain's Seat (K0821 model only)
18" (w) x 18" (d) Full-Size Captain's Seat
20" (w) x 20" (d) Full-Size Captain's Seat

ACCESSORIES

These accessories are sold separately and add functionality and convenience to your power mobility product.



ITEM NUMBERS

TITANLB18CS	18" Low-Back Capt. Seat
TITAN18CS	18" Captain's Seat
TITAN20CS	20" Captain's Seat

WARRANTY

Limited Lifetime on frame, 14 months on electronic controller and motor/brake components, six months on batteries.

SPECIFICATIONS

HCPCS	K0821 (TITANLB18CS) K0823 (TITAN18CS &
	TITAN20CS)
Weight Capacity	
Top Speed	4 mph
Maximum Range	e 15 miles
Turning Radius	25"
Climbing Angle	6°
Ground Clearand	ce 2.5"
Dimensions	40" x 25"
Floor to Seat He	ight 20.5" - 23"
Controller	Dynamic Shark 40A
	24V x 200W x 3900rpm
Batteries (pair)	12V x 35AH (U1)
Charger	3A Offboard
Brakes	
	Electromagnetic
Drive Wheels	10" x 3" Flat-Free
Caster Wheels	6" x 2" Solid
Base Weight	68 lbs.
Battery Weight (pair) 44 lbs.
Seat Weight	39 lbs.
Total Weight	151 lbs.
-	

