



# Jazzy<sup>®</sup>

## 1450

The Jazzy<sup>®</sup> 1450 is a front-wheel drive bariatric power base which has a weight capacity of 600 lbs. and durable construction to make it an exceptionally rugged power chair. Its unique design delivers high-performance operation indoors and out.



**Pride<sup>®</sup>**  
Mobility Products Corp.  
**LIVE YOUR BEST<sup>®</sup>**

## SPECIFICATIONS

<b>Weight capacity</b>	600 lbs.
<b>Maximum speed</b>	Up to 4 mph
<b>Ground clearance</b>	1.5"
<b>Turning radius</b>	26.5"
<b>Overall length</b>	38.25" without front riggings
<b>Overall width</b>	29.25"
<b>Drive wheels</b>	14" solid
<b>Front wheels</b>	4" solid anti-tips
<b>Rear wheels</b>	9" casters
<b>Suspension</b>	Swing Arm Suspension
<b>Drivetrain</b>	Two-motor, front-wheel
<b>Braking system</b>	Intelligent braking (electronic, regenerative disc brakes)
<b>Standard electronics</b>	90 Amp PG VR2
<b>Specialty controls</b>	N/A
<b>Battery charger</b>	Off-board, 8 amp
<b>Per-charge range</b>	Up to 12.14 miles
<b>Battery requirements</b>	Two 12 volt, deep cycle, Group 24, 70 Ah
<b>Battery weight</b>	52 lbs. each
<b>Base weight</b>	195 lbs.
<b>Standard seat weight</b>	37 lbs. (medium-back)
<b>Max. seat dimensions</b>	32" x 28"
<b>Warranty</b>	Lifetime limited warranty on frame; 2-year warranty on electronics; 18-month warranty on drive motors

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

## MODELS

Jazzy 1450-2VHD-SS	(Code: K0826)
Jazzy 1450-2VHD-C	(Code: K0827)
Jazzy 1450-2SPVHD-SS	(Code: K0839)

## FEATURES

- 600 lbs. weight capacity
- Maximum speed up to 4.5 mph
- High-performance, high-efficiency motor package
- Seat dimensions up to 32" x 28"

## OPTIONS

- Synergy® Seat
- High-back seat with solid seat pan
- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder

## COLOR CHOICES



Jazzy Red



Jazzy Blue