

# NEW & IMPROVED!

## CRUISER III

REMOVABLE  
FLIP BACK  
ARMS ALLOW  
FOR EASY  
TRANSFER!



NEW  
TOOL FREE  
ADJUSTABLE  
LENGTH  
FRONT  
RIGGINGS!



Step 1



Step 2



Step 3



ALL THE FEATURES OF THE ORIGINAL  
CRUISER III AND THE NEW BUILT-IN SEAT RAIL  
SYSTEM AND EXTENDABLE UPHOLSTERY!

SEAT RAIL TELESCOPES  
AND UPHOLSTERY  
EXTENDS EASILY  
FROM 16" TO 18"!

**Cruiser III**  
Lightweight, Dual Axle (K0003)

Available  
**1 Pack**

### FEATURES:

- Adjustable seat rail extension and extendable upholstery easily adjust seat depth from 16" to 18" (patent pending)
- Weighs under 36 lbs. (excluding front riggings).
- Carbon steel frame with silver vein finish.
- Removable flip back arms allow for easy transfer.
- New frame style eliminates seat guides and allows for custom back inserts and accessories.
- Nylon upholstery is durable, lightweight, attractive and easy to clean.
- Composite, Mag-style wheels are lightweight and maintenance free.
- 8" front casters are adjustable in 3 positions.
- Padded armrests provide additional comfort.
- Comes with swing-away footrests or elevating legrests with tool free adjustable length riggings.
- Precision sealed wheel bearings in front and rear ensure long lasting performance and reliability.
- Dual axle provides easy transition of seat height to hemi level.
- Comes with push to lock wheel locks.
- Limited Lifetime Warranty.

Due to product improvements and changes, all images in this catalog are subject to change.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.